REQUEST FOR VIRAL SEROLOGY, CULTURE AND MOLECULAR DIAGNOSTICS

New Jersey Department of Health Public Health and Environmental Laboratories 3 Schwarzkopf Drive | Ewing, NJ 08628 https://www.nj.gov/health/phel/

State Lab Specimen ID No.

s.//www.nj.gov/neann/pnei/									Note: Shaded bo	xes are for official use or
Patient Information										
				Sex Male Female Unknown	Date of Birth			CDRSS Number (if applicable)		Outbreak # (if applicable) E-
Patient Address (Street, Apt. #)		State Zip Code Telephone N					umber			
Ethnicity Race Hispanic or Latino White Native Hawaiian or Other Pacific Islander Non-Hispanic or Non-Latino Black or African American American Indian or Alaskan Native Unknown Asian Other								Unknown		
Where (Countries):				hcy Status Hospitalization Status ICU ant Unknown Outpatient Yes egnant Not Applicable Emergency Department No Unknown Unknown Unknown					Admission Date	
Symptom Onset Date:	l Inform	mation (brief history, clinical findings, relevant lab data)								
Relevant Treatment:	Relevant Immunizations: Date:									
Submitting Essilia//		formation				Bogu	ooting D	hygigign	Information	
Submitting Facility/La	Requesting Physician Information Physician Name									
Submitter Address (Street and Suite, City, State, Zip Code)				Physician Address (Street and Suite, City, State, Zip Code)						
Submitter Contact Name				Contact Name						
Contact Telephone No. Contact Fax No.				Contact Tele	Contact Telephone No. Contact Fax No.					
Email Address	Email Addres	Email Address Patient ID No.								
Specimen Information Specimen ID	с	ollection Date			Time			_ 🗆 AI	М 🗌 РМ	NJDOH TEST CODE
Specimen Type Serum Bronchoalveolar Lavage/Wash Stool Plasma (EDTA) Sputum Biopsy/Autopsy CSF Swab (specify) Fixed Tissue Nasal Wash Lesion/Vesicle Aspirate Frozen Tissue Urine Amniotic Fluid Other: NOTE: Zika testing on urine, CSF, and/or amniotic fluid must be accompanied by a serum specimen. Specimen Source Site (i.e., arm, nasopharynx, joint):										
Decimen Source Site (i.e., arm, nasopharynx, joint): est(s) Requested ral Serology Screens M420 Rubella (German Measles) IgG M420 Rubeola (Measles) IgG M430 Rubeola (Measles) IgG M440 Mumps IgG M440 Mumps IgG M440 Mumps IgG M550 Varicella IgG M560 Cytomegalovirus IgG M560 Cytomegalovirus IgG M560 Hepatitis Testing M600 Herpes Group IgG M030 SARS-CoV-2 IgG M600 SARS-CoV-2 IgG M600 Respiratory Viral Surveillance* M650 Respiratory Pathogen Panel M800 MERS CoV M7710 HSV 1 M710 HSV 2 M710 HSV 2 M710 SARS-CoV-2				1 Example 1 Example 1 Example 1 Example 2 Exam	Arboviral Testing 80001 West Nile Virus IgM-Serum 80005 West Nile Virus IgM-CSF 80015 Powassan IgM-Serum 80016 Powassan IgM-CSF 90680 Eastern Equine Encephalitis IgM-Serum 90685 Eastern Equine Encephalitis IgM-CSF 90690 St. Louis Encephalitis IgM-Serum 90691 St. Louis Encephalitis IgM-CSF Zika Testing 80080/82/84 80080/82/84 Zika Trioplex PCR Only 80092/94 Zika Trioplex PCR Only 80090 Follow Up IgM 80050/65 Amnio & Serum PCR * Specimen submission requires prior approval					DATE/TIME RECEIVED

See Instructions for SRD-1 Form

REQUEST FOR VIRAL SEROLOGY, CULTURE AND MOLECULAR DIAGNOSTICS INSTRUCTIONS FOR COMPLETING THE SRD-1 FORM

Please direct clinical or epidemiological questions to the NJ Communicable Disease Service and specimen submission or laboratory questions to the NJ Public Health Laboratory

For All Requests:

- Complete a separate form for <u>each</u> specimen.
- Provide all information requested on the form.
 - Specimens may be rejected and testing will be delayed if information is missing, incomplete, or inaccurate.
 - Please include additional patient information as warranted in the "Pertinent Clinical Information" box on the form, e.g., brief history, clinical findings, relevant lab data, etc.
- Should you need additional space, please submit a separate sheet of paper that includes patient name and date of birth.
- Specimens should be labeled with two (2) identifiers, i.e., Name, date of birth or other unique patient ID.
- Specimen and SRD-1 identifiers must match EXACTLY.
- Ship refrigerated specimens on cold packs and frozen specimens on frozen packs or dry ice for next day delivery.

For Zika Testing:

- For newborn specimen requirements contact the Communicable Disease Service at 609-826-5964.
- For all others: 3 ml of serum and 3 ml of urine are required.

Zika specimens collected and shipped Monday—Thursday

- Refrigerate serum and urine at 2-8°C.
- Ship on frozen cold packs within 24 hours of collection.
- Follow IATA Packaging Instructions for Category B infectious substances.
- If shipment cannot occur within 24 hours, freeze at

⁻15 to -25°C.

Zika specimens collected Friday through Sunday

- Freeze serum and urine at -15 to -25°C.
- Ship Monday on dry ice.
- Follow IATA Packaging instructions for Category B infectious substances and dry ice.
- For more information on PHEL Zika Testing go to: http://www.nj.gov/health/phel/

HSV Molecular Testing

- Collect specimens aseptically as soon as possible after symptom onset or at autopsy.
- Refrigerate specimens immediately and deliver to PHEL as soon as possible.

For Viral Serology:

- Collect acute specimen via venipuncture into appropriate tube (red top, serum separator) within 7 days of onset.
- Convalescent samples should be drawn similarly 14 to 21 days after the acute sample. Allow blood to clot before centrifuging to separate the serum.
- Store sera at 2°C to 8°C (35.6 °F to 46.4°F) for up to seven days until they are shipped to the lab.
- If specimen will not be delivered to the lab within 7 days, freeze serum samples at -20°C (-4°F).

For Molecular/PCR:

Norovirus Testing

- Collect whole stool or Cary Blair medium as soon as possible.
- If specimens will be delivered to the lab within 72 hours, refrigerate at 4°C (39°F).
- If specimens will not be delivered to the lab within 72 hours, freeze at -20°C (-4°F).

<u>Respiratory Specimens</u> (Respiratory Virus Surveillance, Respiratory Pathogen Panel, MERS CoV, Influenza and SARS-CoV-2 RT-PCR Testing):

- Collect respiratory specimens aseptically as soon as possible after symptom onset or at autopsy.
- Refrigerate samples immediately and deliver to the PHEL as soon as possible.
- If delivery will be delayed, specimens should be frozen at -70°C.
- Respiratory Virus Surveillance specimens should be submitted from certified enhanced surveillance providers only.
- All providers may submit samples for influenza testing for subtype identification, or if a sample has been determined inconclusive or unsubtypeable.
- All providers may submit samples for Respiratory Pathogen Panel to aid in the identification of noninfluenza respiratory viruses
- Suspect MERS specimens should be submitted only after consultation with CDS

Package and Ship all specimens according to USDOT Hazardous Materials Transport regulations 49 CFR parts 171-178. Package and ship all Category A infectious substances following USDOT 49 CFR 173.196. Package and ship all Category B infectious substances following the USDOT 49 CFR 173.199. Overnight deliveries are accepted Tuesday through Friday 8:30 am to 3:30 pm; except state holidays. Ground deliveries are accepted Monday through Friday 8:30 am to 3:30 pm; except state holidays. For further questions regarding packaging and shipping please contact PHEL Receiving at (609) 530-8387. Please find link to the state courier service here: https://www.nj.gov/health/phel/.

Ground deliveries should be made to:

New Jersey Public Health and Environmental Laboratories Specimen Receiving Unit 3 Schwarzkopf Drive Ewing, NJ 08628